0318

HLA-B*27

Detection of the HLA-B*27 allele with Real Time PCR

REV ST.RT-53-ENITA.1

REF RT-53 - 48 tests

Instructions For Use

INTRODUCTION AND PURPOSE OF USE

The kit HLA-B*27 is a qualitative test that allows the determination, by means of Real Time PCR, of the presence of HLA-B*27 allele, associated with ankylosing spondylitis (90% of cases).

The Procedure allows the detection of the DNA target by means a genomic amplification reaction. The analysis of the results is made using a Real Time PCR analyzer (thermal cycler integrated with a system for fluorescence detection and a dedicated software).

CONTENT

The kit contains reagents enough to perform 48 amplification tests:

	Quantity	Descrizione
R1	3 x 220 µl	Amplification mMix dNTPs, Tris-HCl, KCl, MgCl ₂ , Taq Polymerase, Nuclease-free water, <i>AmpErase</i> Uracil N-Glycosylase (<i>UNG</i>), ROX (pink Cap)
R2	3 x 130 µl	HLA-B*27 probes Mix Upstream primer 1, Downstream primer 2, downstream primer, Target probes (FAM for HLA-B*27 and VIC internal positive control) Nuclease-free water (white Cap)
R3	3 x 35 µl	HLA-B*27 positive control Cloned DNA corresponding to the HLA-B*27 allele
R4	3 x 35 µl	HLA-B*27 negative control Cloned DNA corresponding to the β-globin gene
R5	1 x 30 µl	Negative Control

MATERIALS AND STRUMENTATION REQUIRED BUT NOT SUPPLIED

Disposable latex powder-free gloves or similar material;
Bench microcentrifuge (12,000 - 14,000 rpm);
Micropipettes and Sterile tips with aerosol filter;
Vortex;
Plastic materials (microplate and optical adhesive cover);
Dry block shaker for 1.5ml conical tubes
Magnetic rack for 1.5ml conical tubes
EZ1 ADV XL DSP DNA Blood Card (ref. 9018702)

Reagents

The **HLA-B*27** kit was developed and validated to be used with the following extraction method:

Manual Extraction

Ref. 51304/51306

QIAmp DNA mini kit. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 50/250 samples. (QIAGEN)

Automatic extraction

Ref. 62124

EZ1 DSP DNA Blood kit. The kit allows the automatic DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 48 samples. (QIAGEN)

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10629800 - *VERSANT® Sample Preparation 1.2 Reagents kit box 1*. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

10629801 - *VERSANT® Sample Preparation 1.2 Reagents kit box 2*. The kit allows the manual DNA extraction from Human samples. The kit contains reagents enough to perform the DNA extraction for 96 samples.

ELITE InGenius®SYSTEM

ELITE InGenius® (ELITechGroup S.p.A., ref. INT030)

For automatic sample analysis with the instrument «ELITE InGenius®» (ELITechGroup S.p.A., ref. INT030) the following generic products are required:

- The extraction cartridges «ELITE InGenius® SP 200» (ELITechGroup S.p.A., ref. INT032SP200)
- The consumables for extraction of nucleic acids from biological samples «ELITE InGenius® SP 200 Consumable Set» (ELITechGroup S.p.A, ref. INT032CS),
- The box for tips waste «ELITE InGenius® Waste Box» (ELITechGroup S.p.A, ref. F2102-000),
- The cassettes for amplification of nucleic acids «ELITE InGenius® PCR Cassette» (ELITechGroup S.p.A, ref. INT035PCR),
- The filter tips for the single nozzle pipettor «300 µL Filter Tips Axygen» (Axygen BioScience Inc., CA, USA, ref. TF-350-L-R-S).

The equipment should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure an optimal performance.

Instruments

The **HLA-B*27** kit was developed and validated to be used with the following instruments:

Extraction System

Ref. 9001492. *EZ1 Advanced XL*.

Robotic Workstation for the automatic purification of the nucleic acids until 14 samples simultaneously (QIAGEN)

Real Time PCR

- 7500 Fast from Lifetechnologies
- StepOne plus from Lifetechnologies
- VERSANT kPCR AD from Siemens or Stratagene MX3005P/MX3000P
- Rotor Gene Q MDx from QIAGEN
- CFX96 Real Time PCR System from Bio-Rad
- LightCycler 480 from Roche

Please ensure that the instruments have been installed, calibrated, checked and maintained according to the manufacturers' instructions and recommendations.

SAMPLES

The **HLA-B*27** system must be used with extracted DNA from the following biological samples: **whole Blood EDTA**. Collected samples must be shipped and stored at +2 - +8°C and used within 3 days from the collected data.

Store the sample at -20°C if it is used after 3 days..

PRECAUTIONS FOR USE

- This kit is for in vitro diagnostics (IVD), for professional use only and not for in vivo use.
- Carefully read this Instruction For Use before using the kit.
- If required, Clonit offers the necessary technical support for the correct use of the kit.
- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood.
- Periodically wipe the working area with 0,5% hypochlorite.
- Wear protective clothing such as laboratory coats and disposable gloves while assaying samples.
- Use powder-free gloves. Do not leave fingerprints on optical caps. Do not write on caps as this may cause an interference with fluorescent detection.
- Avoid any contact between hands and eyes or nose during specimen collection and testing.
- Materials containing or potentially containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood. Waste should be discarded according to local law.
- Never pipette solutions by mouth.
- Use calibrated and regularly checked pipettes and instrumentation only.
- Avoid the air bubbles during the master mix dispensing. Eliminate them before starting amplification.
- Wash hands carefully after handling samples and reagents.
- Do not eat, drink or smoke in the area where specimens and kit reagents are handled.
- Provided reagents are not infectious and hazardous for the health (see Material Safety data Sheet – MSDS).
- The ELITe InGenius® PCR Cassettes must be handled in such a way to reduce as much as possible amplification product diffusion into the environment in order to avoid sample and reagent contamination.

LIMITS OF THE METHOD

The extreme sensitivity of gene amplification may cause false positives due to cross-contamination between samples and/or controls. Therefore, you should:

-physically separate all the products and reagents used for amplification reactions from those used for other reactions, as well as from post-amplification products

- use tips with filters to prevent cross-contamination between samples
- use disposable gloves and change them frequently
- carefully open test tubes to prevent aerosol formation
- close every test tube before opening another one

As with any diagnostic device, the results obtained with this product must be interpreted taking into consideration all the clinical data and other laboratory tests available for the patient

As with any diagnostic device, with this product there is a residual risk of obtaining invalid, false positives or false negatives results.

Patient Drug treatment may interfere with the final result of the molecular biology analysis. The proper functioning of the amplification mix depends on the correct collection, correct transportation, correct storage and correct preparation of a biological sample.

WARNINGS

- Read carefully the instructions for use before using this test.
- Use only extracted DNA from Blood EDTA
- Do not mix reagents from different lots.
- Thaw and carefully mix the reagents of the kit before use.
- The PCR mix has to be freshly prepared every time.

After reconstitution, the amplification master mix must be used in one time (16 reactions). Repeat thawing and freezing of reagents (more than twice) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

- Do not use beyond the expiration date which appears on the package label.
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use.
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request).
- In case of damaged package, contact the technical support before using the kit.
- In case of any serious incident that has occurred in relation to the device, a notice shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

STORAGE AND STABILITY

Store the product **HLA-B*27** at **-10°C/-20°C**.

The **HLA-B*27** kit is shipped on dry ice. The kit components should be frozen. If one or more components are not frozen upon receipt or if the tubes have been compromised during transport, contact Clonit srl for assistance.

An intact and well stored product has a stability of 12 months from the date of production. Do not use beyond the expiration date which appears on the package label.

Repeat thawing and freezing of reagents (more than four times) should be avoided, as this might affect the performance of the assay. The reagents should be frozen in aliquots, if they are to be used intermittently.

The data of the stability tests show that the kit cod RT-53 maintains an optimum amplification efficiency and is stable for 24 Hours when stored at the room temperature and normal light conditions

ANALYTICAL PROCEDURE

Manual Extraction

Ref. 51304/51306 - QIAmp DNA mini kit (QIAGEN).

Follow the instructions inside the kit QIAmp DNA Mini Kit. Elute the sample in 50 µl of buffer AE.

Samples are now ready for amplification or storage at -20°C.

Automatic extraction

Ref. 62124 - EZ1 DSP DNA Blood kit on EZ1 Advanced XL instrument.

Follow the instructions inside the kit EZ1 DSP DNA Blood kit. Start from 200 µl of sample and elute it in 50 µl of buffer AE.

Samples are now ready for amplification or storage at -20°C

Manual extraction (SIEMENS)

Ref. 10629800 - VERSANT® Sample Preparation 1.2 Reagents kit box 1.

Ref. 10629801 - VERSANT® Sample Preparation 1.2 Reagents kit box 2.

Follow the instructions supplied by Siemens and elute it in 70 µl of Elution buffer. Transfer 55 µl of eluted sample to an appropriately size tube.

Samples are now ready for amplification or storage at -20°C

OPERATING PROCEDURE FOR ELITE InGenius®SYSTEM

Before starting the session, referring to the instrument documentation, it is necessary to:

- switch on the ELITE InGenius® and select the login mode "CLOSED",
- verify that the amplification controls (**HLA-B*27 Positive Control**, **HLA-B*27 Negative Control** and **Negative Control**) were run in association with the amplification reagent lot to be used and that they are approved and valid (Status). If there are not amplification controls approved or valid, run them as described in the following paragraphs;
- Choose the type of run, following the instructions on the Graphical User Interface (GUI) for the session setup and using the Assay Protocols provided by ELITechGroup S.p.A. These IVD protocols were specifically validated with the ELITE InGenius® instrument and the cited matrix. The Assay protocol available for sample testing with the **HLA-B*27 Kit** product is described in the table below.

Name	Matrix	Report	Characteristics
Clonit HLA B27_200_200	Whole Blood	Positive Negative	Extraction Input Volume: 200 µL Extraction Elute Volume: 200 µL Internal Control: NO Sonication: NO Dilution Factor: 1 PCR Mix volume: 20 µL Sample PCR input volume: 5 µL

The **HLA-B*27 Kit** product, in association with the ELITE InGenius® system, can be used in order to perform:

- Integrated run (Extract + PCR),
- Amplification run, (PCR only),
- Amplification Positive Control and Negative Control run (PCR only).

All the parameters needed for the session are included in the Assay protocol available on the instrument and are automatically recalled when the Assay protocol is selected.

INTEGRATED RUN (EXTRACT + PCR).

To setup an integrated run, carry out the following steps as per the GUI:

1. Thaw the Amplification Mix and the HLA-B*27 probe mix tube and prepare the PCR master mix by carefully mixing the necessary volume for the number of controls to be tested, as reported in the "Preparation of the PCR mix" section. Prepare the Master Mix in a 2 mL Sarstedt tube (Ref. 72694005).
2. Select "Perform Run" from the "Home" screen.
3. Ensure that the "Extraction Input Volume" is 200 µL and the Extracted Elute Volume is 200 µL.
4. For each Track of interest fill in the "SampleID" (SID) by typing or by scanning the sample barcode.
5. Select the Assay protocol to be used in the "Assay" column
6. Ensure that the "Protocol" displayed is: "Extract + PCR".
7. Select the sample loading position in the "Sample Position" column:
 - if a primary tube is used, select "Primary Tube",
 - if a secondary tube is used, select "Extraction Tube".
8. Click "Next" to continue the setup.
9. Load the master mix prepared in step 1 on the "Inventory Block" selected by following the GUI instruction. Click "Next" button to continue the setup.
10. Load and check the Tip Racks in the "Inventory Area" selected by following the GUI instruction. Click "Next" button to continue the setup.
11. Load the "PCR Cassettes", the "ELITE InGenius® SP 200" extraction cartridges, all the required consumables and the samples to be extracted, following the GUI instruction. Click "Next" to continue the setup.
12. Close the instrument door.
13. Press "Start" to start the run.

After process completion, the ELITE InGenius® system allows users to view, approve, store the results and to print and save the report.

AMPLIFICATION RUN (PCR ONLY)

To set up the amplification run carry out the following steps as per GUI:

1. Thaw the Amplification Mix and the HLA-B*27 probe mix tube and prepare the PCR master mix by carefully mixing the necessary volume for the number of controls to be tested, as reported in the "Preparation of the PCR mix" section. Prepare the Master Mix in a 2 mL Sarstedt tube (Ref. 72694005).
2. Select "Perform Run" from the "Home" screen.
3. Even if no extraction will be carried out, ensure that the Extraction Input Volume is 200 µL and the Extracted Elute Volume is 200 µL.
4. For each Track of interest fill in the "SampleID" (SID) by typing or by scanning the sample barcode.
5. Select the Assay protocol to be used in the "Assay" column
6. Select "PCR Only" in the "Protocol" column.
7. Ensure the sample loading position in the "Sample Position" column is "Elution Tube (bottom row)". Click "Next" to continue the setup.
8. Load the master mix prepared in step 1 on the "Inventory Block" selected by following the GUI instruction. Click "Next" button to continue the setup. Load and check the Tip Racks in the "Inventory Area" selected by following the GUI instruction. Click "Next" to continue the setup.
9. Load the "PCR Cassette" and the extracted Nucleic Acid samples following the GUI instruction. Click "Next" to continue the setup.
10. Close the instrument door.
11. Press "Start" to start the run.

After process completion, the ELITE InGenius® system allows users to view, approve, store the results and to print and save the report.

AMPLIFICATION RUN FOR POSITIVE CONTROL HLA-B*27, NEGATIVE CONTROL HLA-B*27 AND NEGATIVE CONTROL (PCR ONLY)

Before analysis of any sample, it is mandatory to generate and to approve the amplification controls for the amplification reagent lot that will be used in testing:

- as amplification Positive Controls, use the Positive Control HLA-B*27 and the Negative Control HLA-B*27 (provided with this kit) in association with Assay Protocol Clonit HLA B27_PC,
- as amplification Negative Control, use the negative control tube (provided with this kit) in association with Assay Protocol Clonit HLA B27_NC.

To setup the amplification run for amplification controls carry out the following steps as per GUI:

1. Thaw the Amplification Mix and the HLA-B*27 probe mix tube and prepare the PCR master mix by carefully mixing the necessary volume for the number of controls to be tested, as reported in the "Preparation of the PCR mix" section. Prepare the Master Mix in a 2 mL Sarstedt tube (Ref. 72694005).
2. Select "Perform Run" from the "Home" screen.
3. Thaw the Positive Control HLA-B*27 tube, the Negative Control HLA-B*27 tube and the negative control tube for the session.
4. Transfer at least 10 µl of each control to an "Elution tube", provided with the ELITE InGenius® SP 200 Consumable Set.
5. In the Track of interest, select the Assay protocol to be used in the "Assay" column.

- For the Positive Control HLA-B*27 and Negative Control HLA-B*27, select the Assay Protocol Clonit HLA B27_PC in the "Assay" column and fill in the lot number and expiry date
- For the Negative Control, select the Assay Protocol Clonit HLA B27_NC in the "Assay" column and fill in the lot number and expiry date of the molecular biology grade water. Click "Next" to continue the setup.
- Load the master mix prepared in step 1 on the "Inventory Block" selected by following the GUI instruction. Click "Next" button to continue the setup.
- Load and check the Tip Racks in the "Inventory Area" selected by following the GUI instruction. Click "Next" to continue the setup.
- Load the "PCR Cassettes", the Positive Control HLA-B*27 tube, Negative Control HLA-B*27Int and the Negative Control tube following the GUI instruction. Click "Next" to continue the setup.
- Close the instrument door.
- Press "Start" to start the run.

After process completion, the ELITE InGenius® system allows users to view, approve, store the results and to print and save the report.

Review and approval of results

At the end of the run, the "Results Display" screen is automatically shown. In this screen, the sample/control results and the run information are shown. From this screen it is possible to approve the results, print or save the reports ("Sample Report" or "Track Report"). Refer to the ELITE InGenius® instrument user's manual for more details.

The ELITE InGenius® system generates the results with the product **HLA-B*27 Kit** through the following procedure:

- Validation of amplification Positive Control and Negative Control results,
- Validation of sample results,
- Sample result reporting.

The amplification Positive Control HLA-B*27, Negative Control HLA-B*27 and Negative Control results, specific for the amplification reagent lot, will expire **after 15 days**.

The results of Positive Control HLA-B*27, Negative Control HLA-B*27 and Negative Control amplification runs are used by the instrument software to setup the "Control Charts" monitoring the amplification step performances. Refer to the instrument user's manual for more details.

When the Positive Control HLA-B*27 or Negative Control HLA-B*27 or Negative Control result does not meet the acceptance criteria, the "not passed" message is shown on the "Controls" screen and it is not possible to approve it. In this case, the amplification Positive Control HLA-B*27 or Negative Control HLA-B*27 or Negative Control reaction has to be repeated.

When the Positive Control HLA-B*27 and/or Negative Control HLA-B*27 and/or Negative Control is run together with samples to be tested and its result is invalid, the entire session is invalid. In this case, the amplification of all samples must be repeated too.

For each sample, the assay result is automatically interpreted by the system as established by the ELITE InGenius® software algorithm and the Assay protocol parameters. The sample results are stored in the database and can be exported as "Sample Report" and "Track Report".

Sample Result Reporting

The sample results are stored in the database and can be viewed as "Sample Report" and "Track Report".

The "Sample Report" shows the details of a sample run sorted by Sample ID (SID).

The "Track Report" shows the details of a sample run track by track.

The "Sample Report" and "Track Report" can be printed and signed by authorized personnel.

SOFTWARE SETTING

Lifetechnologies 7500 fast/StepOne plus

Turn the instrument and the computer on and open the control software. Click on "**Advance Setup**": by default the software will shows the page "**experiment properties**". Write in the "**experiment name**" the file name, choose the type of instrument (**7500 or 7500fast / StepOne or StepOne Plus**), the type of reaction (**quantitation standard curve**), the type of reagents used (**Taqman®Reagents**) and the analysis reaction time (**Standard ≈ 2 hours to complete a run**).

Open the page named "**page setup**" (sheet **Define Target and Samples**).

In the window **Define Targets** set:

Target	Reporter	Quencer
HLA-B*27 probe:	FAM	NFQ-MGB
IC (β-globin) probe:	VIC	TAMRA

Set the samples' name in the window "**Define Samples**".

In the same page "**plate setup**" select the sheet "**Assign Target and Samples**". On the screen you will see the microplate draft.

Select an area of the plate where the controls will be placed: select wells of the plate and set both targets (HLA-B*27 and β-globin). Select "**Assign target to selected wells**" in the blank, the "**task Standard (S)**" for HLA-B*27 target.

Choose an area in the plate where negative control will be placed: select wells of the plate and set both targets (HLA-B*27 and β-globin). Select "**Assign target to selected wells**" in the blank, the "**task Negative (N)**" for the HLA-B*27 target.

Select an area of the plate where samples will be placed: select the wells and set both targets (HLA-B*27 and β-globine). Link every well to a sample, through the window "**Assign samples to selected wells**".

For each sample, select in the blank "**Assign targets to selected wells**" the "**task UnKnown (U)**" for the HLA-B*27 target.

Set ROX as passive reference, using it as normalizer of detected fluorescence.

Open "**Run Method**" (sheet **Graphic View**) and set the thermal cycling as follows, with the "collect data" in annealing/extension phase:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

In the window "**Reaction volume plate per well**" set a volume of 25 µl.

After having prepared the plate, and correctly inserted in the instrument, press the button "**Start Run**".

Rotor Gene Q MDx

The experiments can be **set using the Quick Start Wizard or the Advanced Wizard**, which appears when the software is started.

Select the wizard "**Advanced**". As a first step, select the model "**Two Step Reaction**" with a double click in the "**New Run**".

In the next window, select the type of rotor installed on the instrument from the list that appears. Check the "**Locking Ring Attached**", check the checkbox and then click "**Next**". Enter the name of the operator and the reaction volume of 25 µl, and then click "**Next**".

In the next window click on "**edit profile**". Set the following thermal cycle:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

Select the annealing/extension from the thermal profile and click on "**Acquiring A to cycling**".

In the next window, select yellow from the **available channels** and add it to **acquiring channel** along with the green channel and click "**OK**". In the next window click on "**OK**" and then click "**Next**".

Click on "Edit Gain" button and set the following values for each channel:

Reporter	Gain
Green	5
Yellow	5

To begin the course, click on the button "Start Run". You can save the model before you begin your run by clicking on "Save Template". After clicking on the button "Start Run" window appears "Save As". The stroke can be saved in the desired position by the user. Once the run started, the window "Edit Samples" allows you to set the name of samples and controls in the positions in which they were loaded on the instrument. Select the locations where they were positioned the controls. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analysed. Select "Positive Control". Enter the concentrations of the controls. Select the location where you placed the Negative Control and name it as Negative Control. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analysed. Select "Negative Controls". Select the location of each sample and enter the name or code of the patient. Clicking on the box next to "Type" correspondent, in the dropdown menu "Samples" you can select the type of sample being analysed. Select "Unknown". At the end of the operation click "OK" in the "edit samples" and wait until the end of the race for the analysis (see "Interpretation of Results").

Versant kPCR AD or Stratagene MX3005P

Turn the instrument on and wait until both green lamps have fixed light, turn on the computer and start the control software. In the main screen will appear the window "New Experiment Options": select "Experiment type": **quantitative PCR (Multiple Standard)**. Turn the lamp on 20 minutes before doing a new experiment. For turning the lamp on, click on the icon of the lamp in the tool bar or select "Lamp On" from the menu "Instruments". Verify the right setting of the gain of the fluorescent reporters: in the menu of settings, choose: "Instrument" and then "Filter set gain setting".

Reporter	Gain
FAM	8
JOE/HEX	4
ROX	1

Click on button "setup" in the tool bar and choose "Plate Setup".

Sign the wells corresponding to positive control. Define the positive control positions in right menu, setting:

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
Pos. Control FAM	FAM/HEX/ROX	ROX	None

Clicking on every single well it will appear the window "well information". Choose the name of the positive control (*HLA-B*27 Positive control*). Identify the wells correspondent to HLA-B*27 Negative control. Define the Negative control positions in the menu by setting:

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
HLA-B*27 Neg Control	FAM/HEX/ROX	ROX	None

Clicking on every single well will appear the window "well information", and you can set HLA-B*27 Neg. Control as the name. Identify the wells correspondent to Negative control. Define the Negative control positions in the menu by setting:

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
NTC	FAM/HEX/ROX	ROX	None

Clicking on every single well will appear the window "well information", and you can set NTC as the name. Sign the wells correspondent to samples. Define the sample positions in the right menu, setting:

Well type:	Collect Fluorescent Data:	Reference Dye:	Replicate Symbol:
Unknown	FAM/HEX/ROX	ROX	None

Clicking on every single well it will appear the window "well information"; it is possible to write in the name or the code of the sample. It's possible to set the name of the dye near the name of the analysed target:

FAM	HEX
HLA-B*27	β -Globin

In the toolbar choose the sheet "Thermal Profile Setup" and set thermal cycle as follows, with the "collect data" in annealing/extension phase:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

After preparing the plate and inserting it in the instrument, press the button "Run", selecting the sheet Thermal profile status and check the accuracy of the thermal profile. Select the box **Turn Lamp Off** at the end of execution. Push the button Start: the software will ask to indicate the name of saved file and will begin the analysis.

CFX 96 Real Time PCR

Turn the instrument and the computer on and start the control software. In the principal screen will appear the window "Startup wizard": select "CFX96" and press "ok". In the next window push "create new" and set the thermal protocol and the reaction volume (25 μ l), with the "collect data" in annealing/extension phase:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95° C 10 min	
35	95° C 15 sec	60° C 1 min

Save the protocol and click the next button. The software will open in default the sheet "plate". Click "create new", select "Fluorophores button" to choose fluorophores (FAM and VIC). Select the locations where they were positioned the controls of known concentration and choose the "Sample Type" **Standards**. Click "Load" check boxes to load fluorophores and Type or select Target Name. Select the location where you placed the Negative Control. Choose the "Sample Type" **NTC**. Click "Load" check boxes to load fluorophores and Type or select Target Name. Select the location of each sample and enter the name or code of the patient. Choose the "Sample Type" **Unknown**. Click "Load" check boxes to load fluorophores and Type or select Target Name.

Save the plate clicking the next button and start the experiment

LightCycler 480

Turn the instrument and the computer on and start the control software. In the principal screen, on "Experiment Creation" select the "Plate type" and push the "New experiment" button. The window "experiment" appears. On the "Run protocol" sheet set: Thermal protocol, Reaction volume (25 μ l) and Detection format (**dual colour hydrolysis probe**), with the "collect data" in annealing/extension phase:

cycles	denaturation	annealing/extension
1	50° C 2 min	
1	95°C 10 min	
35	95° C 15 sec	60° C 1 min

Push the "**Subset editor**" button and in this window, select an area of the plate where controls and samples will be placed.

Push the "**Sample editor**" button. Select the correct workflow (Step1: **Abs quantification**), choose the samples Subset created in the step before and insert the name for each well. Choose the correct type for each well: Positive CTR, Negative CTR or Unknown sample. Push again the "**experiment button**", insert the plate in the instrument and push "**start run**"

PREPARATION OF THE REACTIONS:

Thaw a tube of **Amplification mMix**;

Thaw tube of **HLA-B*27 probe mix**

Mix carefully through vortex **210µl** of **Amplification mMix** with **126µl** of **HLA-B*27 probe mix**

The mix is enough for 16 amplification reactions: **2 positive controls**, **1 negative control** and **13 samples**.

For a total number of tests other than 16, prepare a mixture for $n + 1$ tests, following the volumes shown in the table:

Reagent	Description	Volume for single test (µl)
R1 -	Amplification mMix	12,5
R2	HLA-B*27 probe Mix	7,5

Distribute, in the amplification plate, **20µl** of just reconstituted mix in chosen positions and already setted on the instrument software.

Distribute, in the negative control position, **5µl** of solution taken by **the negative control** vial.

Distribute, in chosen position for each sample, **5µl** of corresponding **sample**.

Distribute, in chosen positions for the positive controls, **5µl** of **HLA-B*27 positive control** and **5µl** of **HLA-B*27 negative control**.

Seal up accurately the plate using an optical adhesive film and verify that there aren't air bubbles in the mix, to avoid the amplification interferences. For the Rotor-Gene Q MDx, seal each tube with the appropriate cap. The air bubbles presence is not influent; the rotor centrifugal force will allow automatic deletion. Transfer the plate in the instrument and push the button "**Start Run**".

QUANTITATIVE ANALYSIS

Lifetechnologies 7500 Fast/StepOne Plus.

At the end of the PCR run, the software automatically opens the "**Analysis**" window in the "**Amplification plot**" sheet on the menu on the left. Select the wells corresponding to the positive control, negative control and samples for analysis.

Select in the "**Option**" window inside the "**Target**" pop-up menu the **HLA-B*27 target**. Check the correct setting of the threshold.

Select in the "**Option**" window inside the "**Target**" pop-up menu the **IC Control target**. Check the correct setting of the threshold.

The analysis of the results is made selecting from the menu in the left the page "**Analysis**". From the page "**Amplification Plot**" verify the amplification plot for every single sample. Opening the sheet "**view well table**" in the right side of the software it is possible to verify the data obtained from experiments: Threshold Cycles, emitted fluorescences etc...

Clicking from the menu file and selecting the box export, the window "**export properties**" will open. Indicate the file name, select the position to save it (**Browse**) and click on button "**Start export**". In this way the software will permit to save a excel file with all the data corresponding to selected experiment.

Rotor-Gene Q MDx

At the end of the PCR run open the "**Analysis**" window. Select the "**Quantification**" sheet and click on "**cycling A (green)**".

Select from the menu "**Dynamic Tube**" and subsequently "**Slope correct**". Check the correct setting of the threshold in the space provided "**CT calculation – Threshold**". Open the "**Analysis**" window. Select the "**Quantification**" sheet and click on "**cycling A (yellow)**". Select from the menu "**Dynamic Tube**" and subsequently "**Slope correct**". Check the correct setting of the threshold in the space provided "**CT calculation – Threshold**". Also in this case, you can print a report of the analysis by clicking on the "**Report**" window and selecting the file in the first **Quantification cycling A (green)** and then the file **cycling A (yellow)**.

Versant kPCR AD or Stratagene MX3005P

Click on button "**Analysis**" in the toolbar. The software will open in default the sheet "**Analysis Term Setting**". Activate the buttons FAM and HEX in the lower part of the screen and select testing samples.

Click on sheet "**results**"; the software will open in default the page "**Amplification plot**". Check the correct setting of the threshold in the specific window "**Threshold fluorescence**", in the menu on the right of the screen. Selecting the box **Text report** from menu "**Area to Analyse**" in the right side of the screen it's possible to verify the data obtained from the experiments: Threshold Cycles, emitted Fluorescences etc.. From the window **Text Report** it's possible to export the results obtained clicking **file, export** on main menu.

CFX96 Real Time PCR System

At the end of the PCR, select the "**quantitation**" sheet. On the top of the screen, select "**settings**" from the menu and choose "**Baseline Threshold...**". You can export the report pushing the paper block figure on the top of the screen.

LightCycler 480

When the run is completed select analysis and choose the correct kind of analysis you want: "**Abs Quant/Fit Points**". Choose the samples subset you want to analyse. Select the "**NoiseBand**" sheet, under the plot you can choose "**NoiseBand (Fluoresc.)**"; and move the line of the NoiseBand on the plot with the mouse of your PC. Repeat this action for each fluorophore using the "**Filter comb**" button. Clicking the sheet "**Analysis**" you can set the threshold choosing the option "Threshold(manual)".

After setting parameters push the "**Calculate**" button. Repeat this action for each fluorophore.

INTERPRETATION OF RESULTS

In the Real Time PCR reaction the Ct values of specific probe for HLA-B*27 are used for detect the presence of the Target in analysis.

Fluorescence increase of the specific probe for HLA-B*27 (FAM) indicates the positivity of the sample for the target in exam.

As with any diagnostic device, the results obtained with this product must be interpreted taking in consideration all the clinical data and other laboratory tests done on the patient.

The use of positive and negative control in each amplification session allows to verify the correct functioning of the amplification mix and the absence of any contamination.

The instrument software is able to analyse the fluorescences that are emitted by the specific probe for HLA-B*27 (FAM) and by the specific probe for the positive internal control (β -globin VIC). A proper functioning of the amplification mix can be verified analysing these parameters:

Parameters	Ref.
Ctrl HLA-B*27 positive (FAM)	Ct \leq 27
Ctrl HLA-B*27 negative (VIC)	Ct \leq 28

If the amplification reaction of each controls produces a Ct > 27 for the CTRL HLA-B*27 positive, and Ct > 28 for the CTRL HLA-B*27 negative, the session can't be considered valid and so it must be cancelled.

Be sure that there isn't any specific fluorescence increasing for examining target in negative control (FAM).

In the amplification reaction of each sample, the Ct values for the internal control (β -globin) specific probe are used for validating the analysis session beginning from extraction process until detection stage.

A good extraction performances presents internal control (β -globin) threshold cycle between 22 and 25.

Be sure that emitted fluorescence from internal control amplification has not a Ct >28 or undetermined. If a sample presents an undetermined HLA-B*27 and internal control Ct >28 means that some problems happened in the extraction stage or in the amplification stage; therefore the sample could be a false negative. **Repeat the sample.**

Detector FAM	Detector VIC/JOE	Assay	Sample
Ct undetermined	Ct > 28 o undetermined	Not valid	Repeat
Ct undetermined	Ct < 28	Valid	Negative
Ct positive	Ct < 28	Valid	Positive
Ct low	Ct > 28 undetermined	Valid	High Positive

PERFORMANCES

Analytical sensitivity:

The highest dilution (titre) to which a positive sample can be diluted without the system losing the ability to detect it as positive with a rate of $\geq 95\%$ is considered analytical sensitivity.

The HLA-B*27 test identifies a congenital target, totally present or absent in the patient's organism. For this reason the analytical sensitivity of the test is not required and the functionality of the system is validated with the correct identification of the internal control (Ct <28).

Clinical sensitivity:

For the purposes of this evaluation is considered as clinical sensitivity the skill of determining true positives in the totality of positive screened samples. The analysis is made on HLA-B*27 positive samples and the test is performed following the method advices. Positive samples are confirmed with another available method.

Positive Samples	Founded Positive Samples	C.sensitivity
69	69	100%

Obtained results show a clinical sensitivity of 100%.

Diagnostic Specificity:

For the purposes of this evaluation is considered as diagnostic specificity the skill of the method of determining real negative samples. The diagnostic specificity of the system is valued analysing human genomic samples tested and confirmed as negative with another available system.

Negative Samples	Founded Negative Samples	D. specificity
112	112	100%

Diagnostic specificity is 100% for material extracted from EDTA blood.

Analytical Specificity:

Test's specificity is guaranteed by the use of specific primers for determining HLA-B*27.

The alignment of the choose regions for specific primers' hybridization for HLA-B*27 demonstrated: their conservation, the absence of significative mutations and the complete specificity for the analysed target.

Samples recognized as positive for a given genotype must be recognized as such by the amplification system described.

Reproducibility and Repeatability:

The reproducibility and repeatability of the system were valued the synthetic positive control of the kit (HLA-B*27 positive control Cloned DNA corresponding to the HLA-B*27 allele and β -globin gene - HLA-B*27 negative control Cloned DNA corresponding to the β -globin gene) quantified by spectrophotometric analysis and a negative control (negative DNA). For each session, 5 replicates were made in 3 different sessions, performed by different technicians, using 3 different lots of product.

Instrument Average Positive CTRL

Positive Control HLA-B*27	Lot	N° replicates	Med. Reveal. Ct. (FAM)	Med. Reveal. Ct. (VIC)	Inaccuracy % (FAM)	Inaccuracy % (VIC)
	L.62	15	21.34	21.34	0.77%	2.38%
	L.63	15	21.32	21.4	1.75%	2.36%
	L.64	15	21.76	21.63	1.97%	1.02%

The medium inaccuracy Positive CTRL HLA-B*27 % is 1,5% and Positive CTRL Beta globine % is 1,92%

Instrument Average Negative CTRL

Negative Control HLA-B*27	Lot	N° replicates	Med. Reveal. Ct.	Inaccuracy %
	L.62	15	21.73	2.37%
	L.63	15	21.67	1.65%
	L.64	15	22.02	0.38%

The medium inaccuracy Negative CTRL Beta globina % is 1,5%.

Reproducibility and Repeatability FOR ELITE InGenius®SYSTEM

To evaluate the Reproducibility and Repeatability of the HLA-B*27 Kit on the ELITE InGenius integrated instrument and to have at least as same performances obtained with other validated instruments.

The ELITE InGenius® system is a fully automated instrument that integrates operations of extraction and purification of nucleic acids, amplification and detection of the target sequence by Real Time polymerase chain reaction (RT-PCR) and interpretation of the results. A dedicated Assay Protocol for the use of the HLA-B*27 Kit on this platform and interpretation of the results were developed in collaboration with ELITechGroup Spa.

To evaluate the reproducibility and repeatability of the test, two contrived samples plus a negative control were tested in 5 replicates using 3 different lots of the product. The contrived samples were comprised of:

- the synthetic positive control R3 of the kit (HLA-B*27 positive control Cloned DNA) corresponding to the HLA-B*27 allele
- the synthetic positive control R4 of the kit (HLA-B*27 negative control Cloned DNA) corresponding to the β -globin gene

The repeatability and reproducibility was performed using Amplification Run (PCR Only)

Instrument Average Positive CTRL

Positive Control HLA-B*27	Lot	N° replicates	Average Ct (FAM)	Average Ct (VIC)	Inaccuracy % FAM	Inaccuracy % VIC
	L.290921	5	20,99	20,26	3,38%	1,55%
	L.161121	5	20,52	20,25	4,77%	3,98%
	L.191221	5	19,38	20,52	1,44%	0,45%

The average inaccuracy for the Positive Control is 2,55 % for the HLA-B*27 allele and 1,99% for the β -globin gene

Instrument Average Negative CTRL

Negative Control HLA-B*27	Lot	N° replicates	Average Ct (VIC)	Inaccuracy %
	L.290921	5	20,60	4,84 %
	L.161121	5	20,48	2,56 %
	L.191221	5	20,68	0,52 %

The average inaccuracy for the Negative Control is 2,64 %

CROSS-REACTIVITY

To check the cross-reactivity of the assay, we use samples tested as negative for HLA-B*27 that meas to have another HLA allele.
We found the 100% of correlation.

INTERFERENCES:

Verify that in the DNA extracted from the sample there is no contamination from mucoproteins and haemoglobin, to exclude possible inhibition of PCR reaction. The interference due to contaminants can be detected through a spectrophotometric analysis, verifying the ratio between the absorbance readings at 260 nm (maximum absorption of Nucleic Acids) and 280 nm (maximum absorption of Proteins). A pure DNA should have a ratio of approximately 1.8.

QUALITY CONTROL

It is recommended to include in each analytical run, as quality control of every extraction, amplification and detection step, an already tested negative and positive sample, or a reference material with known concentration.
Each batch of HLA-B*27 was produced according to predetermined specifications in order to ensure compliance with the ISO EN 13485 certified quality management system.

BIBLIOGRAPHY

Sun Young Cho, M.D., Kwang Gil Lee, M.T., Su Yon Park, M.D., and Hee Joo Lee, M.D. **Utility of In-House PCR for HLA-B27 Typing: Comparison of Concordance Rate between PCR Kit and In-House PCR.** Korean J Lab Med 2008;28:239-43

O. Nathalang, S. Tantimavanich, K. Nillakupt, P. Arnutti & C. Jaruchaimontree. **HLA-B27 testing in Thai patients using the PCR-SSP technique.** Tissue Antigens ISSN 0001-2815

A. Cauli, G. Dessole, M. T. Fiorillo1, A. Vacca, A. Mameli, P. Bitti2, G. Passiu, R. Sorrentino1 and A. Mathieu. **Increased level of HLA-B27 expression in ankylosing spondylitis patients compared with healthy HLA-B27-positive subjects: a possible further susceptibility factor for the development of disease.** Rheumatology 2002;41:1375–1379

TECHNICAL ASSISTANCE

For any question and support please contact our Technical support:

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CLONIT S.r.l.