

Technical Datasheet

Date: 18.02.20

Page 1 of 2

EN ISO 13485 certified		
Article - No:	02-572-8001	Revision: 8
Remarks		single use only

Specification

Product	Centrifuge tube PP with screw cap PE, 50ml, Ø30x115 mm, conical bottom, tube transparent, cap green, graduated, with writing area, np pcr ready, sterile R, max. RCF 20.000g, CE/IVD
Classification	CE/IVD
Dimensions	height: 115,00 mm (+/- 0,2 mm) diameter: 30,00 mm (+/- 0,2 mm)
Colour	tube transp., cap green
Graduation	grad., writ. area
Nominal Capacity	50ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	cold resistant - 80 °C
Autoclavability	Autocl. up to 121°C (no cap)
Centrifugation	max. RCF 20.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004 free of detectable phthalate
Purity Grade	free of PCR inhibitors free of human DNA und RNA free of DNase and RNase free of pyrogens free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	48 months at controlled storage conditions
Packaging	pcs per bag: 25 pcs per carton: 250

Technical Datasheet

Date: 18.02.20

Page 2 of 2

	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	