

Technical Datasheet

Date: 22.02.17

Page 1 of 1



EN ISO 13485 certified

		EN ISO certified
Article - No:	09-400-5061	Revision: 1
Remarks		single use only

Specification

Product	Inoculating loop/needle PS, 1µl, rigid, Ø4x200 mm, 20pcs/peel-pack, dark green, np pcr ready, sterile R
Dimensions	diameter loop: 0,90 mm (+/- 0,1 mm)
	diameter handle: 4,00 mm (+/- 0,1 mm)
	length: 199,00 mm (+/- 0,2 mm)
Colour	dark green
Nominal Capacity	1µl
Material Characteristics	polystyrene (PS)
Temperature range	-50°C - 55°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 1.000
	pcs per carton: 4.000
	pcs per pallet: 160.000
Minimum order qty.	4.000 pcs
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

