

Technical Datasheet

Date: 18.02.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	12-023-0155	Revision: 10
Remarks		single use only

Specification

Product	Transfer pipette PE, 1ml, with micro tip, length 150 mm, transparent
Dimensions	diameter (stem): 4,50 mm (+/- 0,5 mm)
	diameter (bulb): 13,00 mm (+/- 0,5 mm)
	length (stem): 104,50 mm (+/- 0,5 mm)
	length (bulb): 45,50 mm (+/- 0,5 mm)
	total length: 150,00 mm (+/- 10 mm)
Colour	transparent
Nominal Capacity	1ml
Material Characteristics	low density polyethylene (LDPE)
Temperature range	-50°C - 55°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per carton: 500
	pcs per pallet: 90.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	