

Technical Datasheet

Date: 18.02.20

Page 1 of 2

EN ISO 13485 certified		
Article - No:	12-113-8100	Revision: 3
Remarks	Graduation at 0,5ml / 1ml / 1,5ml/ 2ml / 2,5ml / 3ml (for guidance only)	single use only

Specification

Product	Transfer pipette PE, 3ml, single packed, length 150 mm, transparent, grad., sterile R
Dimensions	diameter (stem): 8,00 mm (+/- 0,5 mm)
	outer diameter (orifice): 3,00 mm (+/- 0,2 mm)
	diameter (bulb): 13,00 mm (+/- 0,5 mm)
	length (stem): 109,00 mm (+/- 0,5 mm)
	length (bulb): 41,00 mm (+/- 0,5 mm)
	total length: 150,00 mm (+/- 10 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	3ml
Material Characteristics	low density polyethylene (LDPE)
Temperature range	-50°C - 55°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per cs: 1.500
	pcs per carton: 1.500
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 18.02.20

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

00-000-0001/02.2020

nerbe plus GmbH & Co. KG
Porschestr. 25
D-21423 Winsen/Luhe
Tel: +49 (0)4171 7849 0
Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de
www.nerbe-plus.de
HRA Lüneburg 203240
USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:
NERBE LAB GmbH, Winsen/Luhe
GF: Heike Nerbe
HRB Lüneburg 207694

Deutsche Bank AG
Kto. 0204 41200
IBAN: DE41 2407 0075 0020 4412 00
SWIFT / BIC DEUTDE2H240

Hamburger Sparkasse
Kto. 1387 140831
IBAN: DE56 2005 0550 1387 1408 31
SWIFT / BIC HASPDEHHXXX

