

Technical Datasheet

Date: 18.02.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	12-063-0155	Revision: 5
Remarks	Graduation at 1ml / 2ml / 3ml/ 4ml / 5ml (for guidance only)	single use only

Specification

Product	Transfer pipette PE, 5ml, length 210 mm, transparent, grad.
Dimensions	diameter (stem): 7,60 mm (+/- 0,5 mm)
	diameter (bulb): 16,00 mm (+/- 0,5 mm)
	length (stem): 155,00 mm (+/- 0,5 mm)
	length (bulb): 54,00 mm (+/- 0,5 mm)
	total length: 210,00 mm (+/- 10 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	5ml
Material Characteristics	low density polyethylene (LDPE)
Temperature range	-50°C - 55°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 5.000
	pcs per pallet: 40.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	