

Technical Datasheet

Date: 18.02.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	04-032-0550	Revision: 3
Remarks	Fitting cap: 04-042-0500 and 04-043-0500	single use only

Specification

Product	PCR microcentrifuge tube PP, 0,1ml, without cap, low profile, 8-strips, transparent, np pcr ready max. RCF 20.000g, autocl. up to 121°C, CE/IVD
Classification	CE/IVD
Colour	transparent
Nominal Capacity	0,1ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 20.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per inner cs: 125
	pcs per carton: 1.250
	pcs per pallet: 56.250
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	