

# Technical Datasheet

Date: 18.02.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	04-212-3500	Revision: 8
Remarks	piercable cap	single use only

## Specification

Product	Microcentrifuge tube PP, premium surface, 1,5ml, attached cap, heat-resistant, transparent, grad., writing area, np pcr ready, max. RCF 20.000g, autocl. up to 121°C, CE/IVD
Classification	CE/IVD
Colour	transparent
Graduation	grad., writ. area
Nominal Capacity	1,5ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 20.000g
Product Characteristics	biologically inert
	physiologically harmless
	constantly usable at temperature of 80 °C up to 100 °C
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Packaging	pcs per bag: 200
	pcs per carton: 5.000
	pcs per pallet: 100.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	