

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

|                      |                    |                        |
|----------------------|--------------------|------------------------|
|                      |                    | EN ISO certified       |
| <b>Article - No:</b> | <b>07-642-0000</b> | <b>Revision: 1</b>     |
| <b>Remarks</b>       |                    | <b>single use only</b> |

## Specification

|                                 |  |
|---------------------------------|--|
| <b>Product</b>                  | Filter tip PP, premium surface, 0-100µl, transparent, moulded rings, np premium bio, sterile R, CE/IVD |
| <b>Classification</b>           | CE/IVD   |
| <b>Dimensions</b>               | length: 59,5 mm (+/- 0,2 mm)   |
|                                 | inner diameter: 0,51 mm (+/- 0,2 mm)   |
| <b>Colour</b>                   | transparent  |
| <b>Graduation</b>               | moulded rings  |
| <b>Nominal Capacity</b>         | 0-100µl  |
| <b>Material Characteristics</b> | polypropylene (PP), premium surface  |
|                                 | filter: high density polyethylene (PE)   |
| <b>Product Characteristics</b>  | biologically inert   |
|                                 | physiologically harmless   |
|                                 | no heavy metals according to EC/94/62  |
| <b>Purity Grade</b>             | free of PCR inhibitors   |
|                                 | free of human DNA  |
|                                 | free of RNase and DNase  |
|                                 | free of ATP  |
|                                 | free of pyrogens and endotoxins  |
| <b>Sterilisation Process</b>    | radiation sterilization SAL 10-6   |
| <b>Shelf - Life</b>             | 60 months at controlled storage conditions   |
| <b>Packaging</b>                | pcs per bag: 1.000   |
|                                 | pcs per carton: 10.000   |
|                                 | pcs per pallet: 280.000  |
| <b>Minimum order qty.</b>       | 1.000 pcs  |
| <b>Storage</b>                  | in dry condition: 40% - 60% rel. humidity  |
|                                 | bags and cartons originally closed   |
|                                 | storage temperature +10°C to maximal +30°C   |
|                                 | pallets cannot be stacked  |

Approved by nerbe plus quality department