

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	07-693-8300	Revision: 1
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 100-1.000µl, 96 pcs/rack, transparent, moulded rings, np premium bio, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	length: 81,79 mm (+/- 0,2 mm)
	inner diameter: 0,89 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	100-1.000µl
Material Characteristics	polypropylene (PP), premium surface
	filter: high density polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of ATP
	free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 768
	pcs per carton: 3.840
	pcs per pallet: 107.520
Minimum order qty.	768 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked

Approved by nerbe plus quality department