

# Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	07-662-0000	Revision: 1
Remarks		single use only

## Specification

Product	Filter tip PP, premium surface, 0-200µl, transparent, moulded rings, np premium bio, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	length: 59,5 mm (+/- 0,2 mm) inner diameter: 0,51 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0-200µl
Material Characteristics	polypropylene (PP), premium surface filter: high density polyethylene (PE)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of ATP free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000 pcs per carton: 10.000 pcs per pallet: 280.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked

Approved by nerbe plus quality department