

Technical Datasheet

Date: 23.04.18

Page 1 of 2

		EN ISO certified
Article - No:	07-612-0000	Revision: 3
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 0,1-10µl, extra long, Eppendorf type, transparent, moulded rings np premium bio, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	length: 45,72 mm (+/- 0,2 mm) inner diameter: 0,46 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0,1-10µl
Material Characteristics	polypropylene (PP), premium surface filter: high density polyethylene (PE)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of ATP free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 1.000 pcs per carton: 10.000 pcs per pallet: 450.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable

Technical Datasheet

Date: 23.04.18

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

nerbe plus GmbH
Daimlerstrasse 3
21423 Winsen
Germany
Tel: +49 (0)4171 7849 0
Fax: +49 (0)4171 645 45

E-Mail: info@nerbe-plus.de
Web: www.nerbe-plus.de
Amtsgericht Lüneburg
HRB 110004
USt-IdNr. DE 11681 7171

Deutsche Bank AG
BLZ 240 700 75
Kto. 0204 41200
IBAN: DE41 2407 0075 0020 4412 00
SWIFT / BIC DEUTDE2H240

Commerzbank AG
BLZ 200 800 00
Kto. 1926 47000
IBAN DE43 2008 0000 0192 6470 00
SWIFT / BIC DRESDEFF200
Gläubiger IdNr.: DE33ZZZ00000308514

