

Technical Datasheet

Date: 23.04.18

Page 1 of 2

		EN ISO certified
Article - No:	07-602-0000	Revision: 3
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 0,1-10µl, short, Gilson type, transparent, moulded rings, np premium bio, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	inner diameter: 0,46 mm (+/- 0,2 mm)
	length: 31,20 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0,1-10µl
Material Characteristics	polypropylene (PP), premium surface
	filter: high density polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of ATP
	free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 10.000
	pcs per pallet: 450.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Technical Datasheet

Date: 23.04.18

Page 2 of 2

Approved by nerbe plus quality department