

Quanty HEV

Katalogové číslo:	RT-41
Balení	24 reakcí
Skladování	-20°C

Odkaz do webového katalogu:

<https://allgene.cz/obchod/quanty-hev-24-reakci/>



Quanty HEV je souprava pro detekci a kvantifikaci viru hepatitidy E pomocí Real Time PCR.

Systém quanty HEV je kvantitativní test, který umožňuje amplifikaci a kvantifikaci RNA pomocí Real Time PCR, konkrétně oblasti ORF3 viru hepatitidy E. Analýza výsledků se provádí pomocí analyzátoru real time PCR.

Kit quanty HEV je kompatibilní s následujícími přístroji pro qPCR:

- 7500 Fast Lifetechnologies
- Rotor-Gene Q MDx QIAGEN
- Versant kPCR AD: Siemens nebo Stratagene MX3005P/MX3000P
- LightCycler 480 Roche
- Mic 4
- CFX Biorad

BIBLIOGRAFIE

Sally A. Baylis, Saeko Mizusawa, Yoshiaki Okada, Kay-Martin O. Hanschmann. Collaborative Study to Establish a World Health Organization International Standard for Hepatitis E Virus RNA for Nucleic Acid Amplification Technology (NAT)-Based Assays. WHO/BS/2011.2175

Sally A. Baylis at al. Standardization of HEV Nucleic Acid Amplification Technique based assays: an Initial study to evaluate a panel of HEV Strains and investigate laboratori performance. Journal of Clinical microbiology, Apr. 2011, p1234-1239

Baylis SA, Terao E, Hanschrmann KM, Collaborative Sudy to establish the 1st WHO International Reference Panel for Hepatitis E Virus RNA genotypes for Nucleic Acid Amplification technology (NAT) based assays. WHO Report 2015, WHO/BS/2015.2264