

Technical Datasheet

Date: 18.02.20

Page 1 of 2

EN ISO 13485 certified		
Article - No:	04-043-0500	Revision: 5
Remarks	fits: 04-082-0000, 04-082-0350, 04-082-0550, 04-082-0650, 04-082-9050, 04-082-9056 04-083-0096, 04-083-0150, 04-083-0640, 04-083-0650, 04-032-0500, 04-032-0550	single use only

Specification

Product	Cap for PCR microcentrifuge tubes 0,1 ml & 0,2 ml PP, flat, 8-strips, highly transparent, np pcr ready, max. RCF 20.000g, autocl. up to 121°C, CE/IVD
Classification	CE/IVD
Colour	highly transparent
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 20.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors free of human DNA und RNA free of DNase and RNase free of pyrogens free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 125 pcs per carton: 1.250 pcs per pallet: 243.750
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable

Technical Datasheet

Date: 18.02.20

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

00-000-0001/02.2020

nerbe plus GmbH & Co. KG
Porschestra. 25
D-21423 Winsen/Luhe
Tel: +49 (0)4171 7849 0
Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de
www.nerbe-plus.de
HRA Lüneburg 203240
USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:
NERBE LAB GmbH, Winsen/Luhe
GF: Heike Nerbe
HRB Lüneburg 207694

Deutsche Bank AG
Kto. 0204 41200
IBAN: DE41 2407 0075 0020 4412 00
SWIFT / BIC DEUTDE2H240

Hamburger Sparkasse
Kto. 1387 140831
IBAN: DE56 2005 0550 1387 1408 31
SWIFT / BIC HASPDEHHXXX

